

Antibodies Product Datasheet

Product Name:	Mad 4 Antibody
Catalog Number:	C11039
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAD4. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	MXD4
Protein Name:	Max dimerization protein 4
Other Name:	MXD4; BHLHC12; MAD4; Max dimerization protein 4; Max dimerizer 4; Class C basic helix-loop-helix protein 12; bHLHc12; Max-associated protein 4; Max-interacting transcriptional repressor MAD4
Human Swiss Prot No.:	Q14582
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Mad 4 Polyclonal Antibody detects endogenous levels of Mad 4 protein.
Storage and Stability:	-20°C/1 year