

Antibodies Product Datasheet

Product Name:	KIAA0100 Antibody
Catalog Number:	C11042
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KIAA0100. AA range:681-730
Concentration:	1 mg/ml
Gene Name:	KIAA0100
Protein Name:	UPF0378 protein KIAA0100
Other Name:	KIAA0100; BCOX1; UPF0378 protein KIAA0100; Antigen MLAA-22; Breast cancer-overexpressed gene 1 protein
Human Swiss Prot No.:	Q14667
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KIAA0100 Polyclonal Antibody detects endogenous levels of KIAA0100 protein.
Storage and Stability:	-20°C/1 year