

Antibodies Product Datasheet

Product Name:	Mi2- β Antibody
Catalog Number:	C11048
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CHD4. AA range:571-620
Concentration:	1 mg/ml
Gene Name:	CHD4
Protein Name:	Chromodomain-helicase-DNA-binding protein 4
Other Name:	CHD4; Chromodomain-helicase-DNA-binding protein 4; CHD-4; ATP-dependent helicase CHD4; Mi-2 autoantigen 218 kDa protein; Mi2-beta
Human Swiss Prot No.:	Q14839
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Mi2- β Polyclonal Antibody detects endogenous levels of Mi2- β protein.
Storage and Stability:	-20°C/1 year