

Antibodies Product Datasheet

Product Name:	IRF-4 Antibody
Catalog Number:	C11070
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IRF4. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	IRF4
Protein Name:	Interferon regulatory factor 4
Other Name:	IRF4; MUM1; Interferon regulatory factor 4; IRF-4; Lymphocyte-specific interferon regulatory factor; LSIRF; Multiple myeloma oncogene 1; NF-EM5
Human Swiss Prot No.:	Q15306
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IRF-4 Polyclonal Antibody detects endogenous levels of IRF-4 protein.
Storage and Stability:	-20°C/1 year