

Antibodies Product Datasheet

Product Name:	BS69 Antibody
Catalog Number:	C11071
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZMY11. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	ZMYND11
Protein Name:	Zinc finger MYND domain-containing protein 11
Other Name:	ZMYND11; BS69; Zinc finger MYND domain-containing protein 11; Adenovirus 5 E1A-binding protein; Protein BS69
Human Swiss Prot No.:	Q15326
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BS69 Polyclonal Antibody detects endogenous levels of BS69 protein.
Storage and Stability:	-20°C/1 year