

Antibodies Product Datasheet

Product Name:	ASC-1 Antibody
Catalog Number:	C11094
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TRIP4. AA range:81-130
Concentration:	1 mg/ml
Gene Name:	TRIP4
Protein Name:	Activating signal cointegrator 1
Other Name:	TRIP4; Activating signal cointegrator 1; ASC-1; Thyroid receptor-interacting protein 4; TR-interacting protein 4; TRIP-4
Human Swiss Prot No.:	Q15650
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ASC-1 Polyclonal Antibody detects endogenous levels of ASC-1 protein.
Storage and Stability:	-20°C/1 year