

Antibodies Product Datasheet

Product Name:	GPR68 Antibody
Catalog Number:	C11103
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SPR1. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	GPR68
Protein Name:	Ovarian cancer G-protein coupled receptor 1
Other Name:	GPR68; OGR1; Ovarian cancer G-protein coupled receptor 1; OGR-1; G-protein coupled receptor 68; GPR12A; Sphingosylphosphorylcholine receptor
Human Swiss Prot No.:	Q15743
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR68 Polyclonal Antibody detects endogenous levels of GPR68 protein.
Storage and Stability:	-20°C/1 year