

Antibodies Product Datasheet

Product Name:	ATBF1 Antibody
Catalog Number:	C11114
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZFH3. AA range:761-810
Concentration:	1 mg/ml
Gene Name:	ZFH3
Protein Name:	Zinc finger homeobox protein 3
Other Name:	ZFH3; ATBF1; Zinc finger homeobox protein 3; AT motif-binding factor; AT-binding transcription factor 1; Alpha-fetoprotein enhancer-binding protein; Zinc finger homeodomain protein 3; ZFH-3
Human Swiss Prot No.:	Q15911
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATBF1 Polyclonal Antibody detects endogenous levels of ATBF1 protein.
Storage and Stability:	-20°C/1 year