

Antibodies Product Datasheet

Product Name:	B-ATF Antibody
Catalog Number:	C11122
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BATF. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	BATF
Protein Name:	Basic leucine zipper transcriptional factor ATF-like
Other Name:	BATF; Basic leucine zipper transcriptional factor ATF-like; B-cell-activating transcription factor; B-ATF; SF-HT-activated gene 2 protein; SFA-2
Human Swiss Prot No.:	Q16520
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	B-ATF Polyclonal Antibody detects endogenous levels of B-ATF protein.
Storage and Stability:	-20°C/1 year