

Antibodies Product Datasheet

Product Name:	E4BP4 Antibody
Catalog Number:	C11134
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NFIL3. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	NFIL3
Protein Name:	Nuclear factor interleukin-3-regulated protein
Other Name:	NFIL3; E4BP4; IL3BP1; Nuclear factor interleukin-3-regulated protein; E4 promoter-binding protein 4; Interleukin-3 promoter transcriptional activator; Interleukin-3-binding protein 1; Transcriptional activator NF-IL3A
Human Swiss Prot No.:	Q16649
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	E4BP4 Polyclonal Antibody detects endogenous levels of E4BP4 protein.
Storage and Stability:	-20°C/1 year