

Antibodies Product Datasheet

Product Name:	Dok-7 Antibody
Catalog Number:	C11148
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DOK7. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	DOK7
Protein Name:	Protein Dok-7
Other Name:	DOK7; C4orf25; Protein Dok-7; Downstream of tyrosine kinase 7
Human Swiss Prot No.:	Q18PE1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dok-7 Polyclonal Antibody detects endogenous levels of Dok-7 protein.
Storage and Stability:	-20°C/1 year