

Antibodies Product Datasheet

Product Name:	LRRK1 Antibody
Catalog Number:	C11153
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LRRK1. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	LRRK1
Protein Name:	Leucine-rich repeat serine/threonine-protein kinase 1
Other Name:	LRRK1; KIAA1790; Leucine-rich repeat serine/threonine-protein kinase 1
Human Swiss Prot No.:	Q38SD2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LRRK1 Polyclonal Antibody detects endogenous levels of LRRK1 protein.
Storage and Stability:	-20°C/1 year