

Antibodies Product Datasheet

Product Name:	DGK- κ Antibody
Catalog Number:	C11181
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DGKK. AA range:1221-1270
Concentration:	1 mg/ml
Gene Name:	DGKK
Protein Name:	Diacylglycerol kinase kappa
Other Name:	DGKK; Diacylglycerol kinase kappa; DAG kinase kappa; 142 kDa diacylglycerol kinase; Diglyceride kinase kappa; DGK-kappa
Human Swiss Prot No.:	Q5KSL6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DGK- κ Polyclonal Antibody detects endogenous levels of DGK- κ protein.
Storage and Stability:	-20°C/1 year