

Antibodies Product Datasheet

Product Name:	Sideroflexin-4 Antibody
Catalog Number:	C11227
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SFXN4. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	SFXN4
Protein Name:	Sideroflexin-4
Other Name:	SFXN4; BCRM1; Sideroflexin-4; Breast cancer resistance marker 1
Human Swiss Prot No.:	Q6P4A7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Sideroflexin-4 Polyclonal Antibody detects endogenous levels of Sideroflexin-4 protein.
Storage and Stability:	-20°C/1 year