

Antibodies Product Datasheet

Product Name:	FAM111B Antibody
Catalog Number:	C11231
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/200. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human F111B. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	FAM111B
Protein Name:	Protein FAM111B
Other Name:	FAM111B; CANP; Protein FAM111B; Cancer-associated nucleoprotein
Human Swiss Prot No.:	Q6SJ93
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FAM111B Polyclonal Antibody detects endogenous levels of FAM111B protein.
Storage and Stability:	-20°C/1 year