

Antibodies Product Datasheet

Product Name:	Emp Antibody
Catalog Number:	C11253
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAEA. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	MAEA
Protein Name:	Macrophage erythroblast attacher
Other Name:	MAEA; EMP; HLC10; PIG5; Macrophage erythroblast attacher; Cell proliferation-inducing gene 5 protein; Erythroblast macrophage protein; Human lung cancer oncogene 10 protein; HLC-10
Human Swiss Prot No.:	Q7L5Y9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Emp Polyclonal Antibody detects endogenous levels of Emp protein.
Storage and Stability:	-20°C/1 year