

Antibodies Product Datasheet

Product Name:	Mob1A Antibody
Catalog Number:	C11256
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MOL1A. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	MOB1B
Protein Name:	MOB kinase activator 1B
Other Name:	MOB1B; MOB4A; MOBKL1A; MOB kinase activator 1B; Mob1 homolog 1A; Mob1A; Mob1B; Mps one binder kinase activator-like 1A
Human Swiss Prot No.:	Q7L9L4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Mob1A Polyclonal Antibody detects endogenous levels of Mob1A protein.
Storage and Stability:	-20°C/1 year