

Antibodies Product Datasheet

Product Name:	ATPBD3 Antibody
Catalog Number:	C11268
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATPBD3. AA range:291-340
Concentration:	1 mg/ml
Gene Name:	CTU1
Protein Name:	Cytoplasmic tRNA 2-thiolation protein 1
Other Name:	CTU1; ATPBD3; NCS6; Cytoplasmic tRNA 2-thiolation protein 1; ATP-binding domain-containing protein 3; Cancer-associated gene protein; Cytoplasmic tRNA adenylyltransferase 1
Human Swiss Prot No.:	Q7Z7A3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATPBD3 Polyclonal Antibody detects endogenous levels of ATPBD3 protein.
Storage and Stability:	-20°C/1 year