

Antibodies Product Datasheet

Product Name:	Alpha-protein Kinase 2 Antibody
Catalog Number:	C11273
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ALPK2. AA range:1741-1790
Concentration:	1 mg/ml
Gene Name:	ALPK2
Protein Name:	Alpha-protein kinase 2
Other Name:	ALPK2; HAK; Alpha-protein kinase 2; Heart alpha-protein kinase
Human Swiss Prot No.:	Q86TB3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Alpha-protein Kinase 2 Polyclonal Antibody detects endogenous levels of Alpha-protein Kinase 2 protein.
Storage and Stability:	-20°C/1 year