

Antibodies Product Datasheet

Product Name:	MAGI-2 Antibody
Catalog Number:	C11278
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAGI2. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	MAGI2
Protein Name:	Membrane-associated guanylate kinase WW and PDZ domain-containing protein 2
Other Name:	MAGI2; ACVRINP1; AIP1; KIAA0705; Membrane-associated guanylate kinase; WW and PDZ domain-containing protein 2; Atrophin-1-interacting protein 1; AIP-1; Atrophin-1-interacting protein A; Membrane-associated guanylate kinase inverted 2; MAGI-Q86UL8
Human Swiss Prot No.:	Q86UL8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MAGI-2 Polyclonal Antibody detects endogenous levels of MAGI-2 protein.
Storage and Stability:	-20°C/1 year