

Antibodies Product Datasheet

Product Name:	GLK Antibody
Catalog Number:	C11306
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAP4K3. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	MAP4K3
Protein Name:	Mitogen-activated protein kinase kinase kinase 3
Other Name:	MAP4K3; RAB8IPL1; Mitogen-activated protein kinase kinase kinase 3; Germinal center kinase-related protein kinase; GLK; MAPK/ERK kinase kinase kinase 3; MEK kinase kinase 3; MEKKK 3
Human Swiss Prot No.:	Q8IVH8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GLK Polyclonal Antibody detects endogenous levels of GLK protein.
Storage and Stability:	-20°C/1 year