

Antibodies Product Datasheet

Product Name:	NKIAMRE Antibody
Catalog Number:	C11309
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDKL3. AA range:291-340
Concentration:	1 mg/ml
Gene Name:	CDKL3
Protein Name:	Cyclin-dependent kinase-like 3
Other Name:	CDKL3; NKIAMRE; Cyclin-dependent kinase-like 3; Serine/threonine-protein kinase NKIAMRE
Human Swiss Prot No.:	Q8IVW4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NKIAMRE Polyclonal Antibody detects endogenous levels of NKIAMRE protein.
Storage and Stability:	-20°C/1 year