

Antibodies Product Datasheet

Product Name:	TP53INP2 Antibody
Catalog Number:	C11320
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TP53INP2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	TP53INP2
Protein Name:	Tumor protein p53-inducible nuclear protein 2
Other Name:	TP53INP2; C20orf110; DOR; PINH; Tumor protein p53-inducible nuclear protein 2; Diabetes and obesity-regulated gene; p53-inducible protein U; PIG-U
Human Swiss Prot No.:	Q8IXH6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TP53INP2 Polyclonal Antibody detects endogenous levels of TP53INP2 protein.
Storage and Stability:	-20°C/1 year