

Antibodies Product Datasheet

Product Name:	Fucokinase Antibody
Catalog Number:	C11334
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FUK. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	FUK
Protein Name:	L-fucose kinase
Other Name:	FUK; L-fucose kinase; Fucokinase
Human Swiss Prot No.:	Q8N0W3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Fucokinase Polyclonal Antibody detects endogenous levels of Fucokinase protein.
Storage and Stability:	-20°C/1 year