

Antibodies Product Datasheet

Product Name:	DBC-1 Antibody
Catalog Number:	C11338
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KIAA1967. AA range:431-480
Concentration:	1 mg/ml
Gene Name:	KIAA1967
Protein Name:	DBIRD complex subunit KIAA1967
Other Name:	KIAA1967; DBC1; DBIRD complex subunit KIAA1967; Deleted in breast cancer gene 1 protein; DBC-1; DBC.1; p30 DBC
Human Swiss Prot No.:	Q8N163
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DBC-1 Polyclonal Antibody detects endogenous levels of DBC-1 protein.
Storage and Stability:	-20°C/1 year