

Antibodies Product Datasheet

Product Name:	Spi-C Antibody
Catalog Number:	C11346
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SPIC. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	SPIC
Protein Name:	Transcription factor Spi-C
Other Name:	SPIC; Transcription factor Spi-C
Human Swiss Prot No.:	Q8N5J4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Spi-C Polyclonal Antibody detects endogenous levels of Spi-C protein.
Storage and Stability:	-20°C/1 year