

Antibodies Product Datasheet

Product Name:	MLKL Antibody
Catalog Number:	C11361
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from MLKL . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	MLKL
Protein Name:	Mixed lineage kinase domain-like protein
Other Name:	MLKL; Mixed lineage kinase domain-like protein
Human Swiss Prot No.:	Q8NB16
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MLKL Polyclonal Antibody detects endogenous levels of MLKL protein.
Storage and Stability:	-20°C/1 year