

Antibodies Product Datasheet

Product Name:	AGR3 Antibody
Catalog Number:	C11390
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AGR3. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	AGR3
Protein Name:	Anterior gradient protein 3 homolog
Other Name:	AGR3; BCMP11; Anterior gradient protein 3 homolog; AG-3; AG3; hAG-3; Breast cancer membrane protein 11
Human Swiss Prot No.:	Q8TD06
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AGR3 Polyclonal Antibody detects endogenous levels of AGR3 protein.
Storage and Stability:	-20°C/1 year