

Antibodies Product Datasheet

Product Name:	ERK 8 Antibody
Catalog Number:	C11391
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAPK15. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	MAPK15
Protein Name:	Mitogen-activated protein kinase 15
Other Name:	MAPK15; ERK7; ERK8; Mitogen-activated protein kinase 15; MAP kinase 15; MAPK 15; Extracellular signal-regulated kinase 7; ERK-7; Extracellular signal-regulated kinase 8; ERK-8
Human Swiss Prot No.:	Q8TD08
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ERK 8 Polyclonal Antibody detects endogenous levels of ERK 8 protein.
Storage and Stability:	-20°C/1 year