

Antibodies Product Datasheet

Product Name:	Eps8L3 Antibody
Catalog Number:	C11394
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EPS8L3. AA range:401-450
Concentration:	1 mg/ml
Gene Name:	EPS8L3
Protein Name:	Epidermal growth factor receptor kinase substrate 8-like protein 3
Other Name:	EPS8L3; EPS8R3; Epidermal growth factor receptor kinase substrate 8-like protein 3; EPS8-like protein 3; Epidermal growth factor receptor pathway substrate 8-related protein 3; EPS8-related protein 3
Human Swiss Prot No.:	Q8TE67
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Eps8L3 Polyclonal Antibody detects endogenous levels of Eps8L3 protein.
Storage and Stability:	-20°C/1 year