

Antibodies Product Datasheet

Product Name:	Eps8L1 Antibody
Catalog Number:	C11395
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ES8L1. AA range:411-460
Concentration:	1 mg/ml
Gene Name:	EPS8L1
Protein Name:	Epidermal growth factor receptor kinase substrate 8-like protein 1
Other Name:	EPS8L1; DRC3; EPS8R1; PP10566; Epidermal growth factor receptor kinase substrate 8-like protein 1; EPS8-like protein 1; Epidermal growth factor receptor pathway substrate 8-related protein 1; EPS8-related protein 1
Human Swiss Prot No.:	Q8TE68
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Eps8L1 Polyclonal Antibody detects endogenous levels of Eps8L1 protein.
Storage and Stability:	-20°C/1 year