

Antibodies Product Datasheet

Product Name:	DUSP19 Antibody
Catalog Number:	C11401
Reactivity:	Human, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DUSP19. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	DUSP19
Protein Name:	Dual specificity protein phosphatase 19
Other Name:	DUSP19; DUSP17; LMWDSP3; SKRP1; Dual specificity protein phosphatase 19; Dual specificity phosphatase TS-DSP1; Low molecular weight dual specificity phosphatase 3; LMW-DSP3; Protein phosphatase SKRP1; Stress-activated protein kinase pathway
Human Swiss Prot No.:	Q8WTR2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DUSP19 Polyclonal Antibody detects endogenous levels of DUSP19 protein.
Storage and Stability:	-20°C/1 year