

Antibodies Product Datasheet

Product Name:	PRPF31 Antibody
Catalog Number:	C11412
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PRP31. AA range:331-380
Concentration:	1 mg/ml
Gene Name:	PRPF31
Protein Name:	U4/U6 small nuclear ribonucleoprotein Prp31
Other Name:	PRPF31; PRP31; U4/U6 small nuclear ribonucleoprotein Prp31; Pre-mRNA-processing factor 31; Serologically defined breast cancer antigen NY-BR-99; U4/U6 snRNP 61 kDa protein; Protein 61K; hPrp31
Human Swiss Prot No.:	Q8WWY3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PRPF31 Polyclonal Antibody detects endogenous levels of PRPF31 protein.
Storage and Stability:	-20°C/1 year