

Antibodies Product Datasheet

Product Name:	CNK2 Antibody
Catalog Number:	C11415
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CNKR2. AA range:641-690
Concentration:	1 mg/ml
Gene Name:	CNKSR2
Protein Name:	Connector enhancer of kinase suppressor of ras 2
Other Name:	CNKSR2; CNK2; KIAA0902; KSR2; Connector enhancer of kinase suppressor of ras 2; Connector enhancer of KSR 2; CNK homolog protein 2; CNK2
Human Swiss Prot No.:	Q8WXI2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CNK2 Polyclonal Antibody detects endogenous levels of CNK2 protein.
Storage and Stability:	-20°C/1 year