

Antibodies Product Datasheet

Product Name:	KG19 Antibody
Catalog Number:	C11418
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BRI3B. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	BRI3BP
Protein Name:	BRI3-binding protein
Other Name:	BRI3BP; KG19; BRI3-binding protein; I3-binding protein; Cervical cancer 1 proto-oncogene-binding protein KG19; HCCRB-1
Human Swiss Prot No.:	Q8WY22
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KG19 Polyclonal Antibody detects endogenous levels of KG19 protein.
Storage and Stability:	-20°C/1 year