

Antibodies Product Datasheet

Product Name:	AKAP 149 Antibody
Catalog Number:	C11427
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AKAP1. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	AKAP1
Protein Name:	A-kinase anchor protein 1 mitochondrial
Other Name:	AKAP1; AKAP149; PRKA1; A-kinase anchor protein 1; mitochondrial; A-kinase anchor protein 149 kDa; AKAP 149; Dual specificity A-kinase-anchoring protein 1; D-AKAP-1; Protein kinase A-anchoring protein 1; PRKA1; Spermatid A-kinase anchor prot
Human Swiss Prot No.:	Q92667
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AKAP 149 Polyclonal Antibody detects endogenous levels of AKAP 149 protein.
Storage and Stability:	-20°C/1 year