

Antibodies Product Datasheet

Product Name:	AP-2 γ Antibody
Catalog Number:	C11429
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AP2C. AA range:401-450
Concentration:	1 mg/ml
Gene Name:	TFAP2C
Protein Name:	Transcription factor AP-2 gamma
Other Name:	TFAP2C; Transcription factor AP-2 gamma; AP2-gamma; Activating enhancer-binding protein 2 gamma; Transcription factor ERF-1
Human Swiss Prot No.:	Q92754
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AP-2 γ Polyclonal Antibody detects endogenous levels of AP-2 γ protein.
Storage and Stability:	-20°C/1 year