

Antibodies Product Datasheet

Product Name:	ENX-2 Antibody
Catalog Number:	C11435
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EZH1. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	EZH1
Protein Name:	Histone-lysine N-methyltransferase EZH1
Other Name:	EZH1; KIAA0388; Histone-lysine N-methyltransferase EZH1; ENX-2; Enhancer of zeste homolog 1
Human Swiss Prot No.:	Q92800
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ENX-2 Polyclonal Antibody detects endogenous levels of ENX-2 protein.
Storage and Stability:	-20°C/1 year