

Antibodies Product Datasheet

Product Name:	CUG-BP1 Antibody
Catalog Number:	C11440
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CELF-1. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	CELF1
Protein Name:	CUGBP Elav-like family member 1
Other Name:	CELF1; BRUNOL2; CUGBP; CUGBP1; NAB50; CUGBP Elav-like family member 1; CELF-1; 50 kDa nuclear polyadenylated RNA-binding protein; Bruno-like protein 2; CUG triplet repeat RNA-binding protein 1; CUG-BP1; CUG-BP- and ETR-3-like factor 1; Dead
Human Swiss Prot No.:	Q92879
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CUG-BP1 Polyclonal Antibody detects endogenous levels of CUG-BP1 protein.
Storage and Stability:	-20°C/1 year