

Antibodies Product Datasheet

Product Name:	TRAP230 Antibody
Catalog Number:	C11455
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MED12. AA range:611-660
Concentration:	1 mg/ml
Gene Name:	MED12
Protein Name:	Mediator of RNA polymerase II transcription subunit 12
Other Name:	MED12; ARC240; CAGH45; HOPA; KIAA0192; TNRC11; TRAP230; Mediator of RNA polymerase II transcription subunit 12; Activator-recruited cofactor 240 kDa component; ARC240; CAG repeat protein 45; Mediator complex subunit 12; OPA-containing prote
Human Swiss Prot No.:	Q93074
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TRAP230 Polyclonal Antibody detects endogenous levels of TRAP230 protein.
Storage and Stability:	-20°C/1 year