

Antibodies Product Datasheet

Product Name:	CNK1 Antibody
Catalog Number:	C11460
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CNKR1. AA range:491-540
Concentration:	1 mg/ml
Gene Name:	CNKSR1
Protein Name:	Connector enhancer of kinase suppressor of ras 1
Other Name:	CNKSR1; CNK1; Connector enhancer of kinase suppressor of ras 1; Connector enhancer of KSR 1; CNK homolog protein 1; CNK1; hCNK1; Connector enhancer of KSR-like
Human Swiss Prot No.:	Q969H4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CNK1 Polyclonal Antibody detects endogenous levels of CNK1 protein.
Storage and Stability:	-20°C/1 year