

Antibodies Product Datasheet

Product Name:	Jamip2 Antibody
Catalog Number:	C11472
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human JAKMIP2. AA range:761-810
Concentration:	1 mg/ml
Gene Name:	JAKMIP2
Protein Name:	Janus kinase and microtubule-interacting protein 2
Other Name:	JAKMIP2; JAMIP2; KIAA0555; NECC1; Janus kinase and microtubule-interacting protein 2; CTCL tumor antigen HD-CL-04; Neuroendocrine long coiled-coil protein 1
Human Swiss Prot No.:	Q96AA8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Jamip2 Polyclonal Antibody detects endogenous levels of Jamip2 protein.
Storage and Stability:	-20°C/1 year