

Antibodies Product Datasheet

Product Name:	Mob3A/B Antibody
Catalog Number:	C11480
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Mob3A/B . at AA range: 1-80
Concentration:	1 mg/ml
Gene Name:	MMOBKL2A/MMOBKL2B
Protein Name:	Mps one binder kinase activator-like 2A/B
Other Name:	MOB3A; MOBKL2A; MOB kinase activator 3A; MOB-LAK; Mob1 homolog 2A; Mps one binder kinase activator-like 2A; MOB3B; MOBKL2B; MOB kinase activator 3B; Mob1 homolog 2b; Mps one binder kinase activator-like 2B
Human Swiss Prot No.:	Q96BX8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Mob3A/B Polyclonal Antibody detects endogenous levels of Mob3A/B protein.
Storage and Stability:	-20°C/1 year