

Antibodies Product Datasheet

Product Name:	Mic-1 Antibody
Catalog Number:	C11486
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MIC1. AA range:591-640
Concentration:	1 mg/ml
Gene Name:	C18orf8
Protein Name:	Uncharacterized protein C18orf8
Other Name:	C18orf8; MIC1; Uncharacterized protein C18orf8; Colon cancer-associated protein Mic1; Mic-1
Human Swiss Prot No.:	Q96DM3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Mic-1 Polyclonal Antibody detects endogenous levels of Mic-1 protein.
Storage and Stability:	-20°C/1 year