

Antibodies Product Datasheet

Product Name:	MARK4 Antibody
Catalog Number:	C11511
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MARK4. AA range:461-510
Concentration:	1 mg/ml
Gene Name:	MARK4
Protein Name:	MAP/microtubule affinity-regulating kinase 4
Other Name:	MARK4; KIAA1860; MARKL1; MAP/microtubule affinity-regulating kinase 4; MAP/microtubule affinity-regulating kinase-like 1
Human Swiss Prot No.:	Q96L34
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MARK4 Polyclonal Antibody detects endogenous levels of MARK4 protein.
Storage and Stability:	-20°C/1 year