

Antibodies Product Datasheet

Product Name:	LMTK3 Antibody
Catalog Number:	C11528
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LMTK3. AA range:1251-1300
Concentration:	1 mg/ml
Gene Name:	LMTK3
Protein Name:	Serine/threonine-protein kinase LMTK3
Other Name:	LMTK3; KIAA1883; TYKLM3; Serine/threonine-protein kinase LMTK3; Lemur tyrosine kinase 3
Human Swiss Prot No.:	Q96Q04
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LMTK3 Polyclonal Antibody detects endogenous levels of LMTK3 protein.
Storage and Stability:	-20°C/1 year