

Antibodies Product Datasheet

Product Name:	ALS2CR7 Antibody
Catalog Number:	C11529
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AL2S7. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	CDK15
Protein Name:	Cyclin-dependent kinase 15
Other Name:	CDK15; ALS2CR7; PFTK2; Cyclin-dependent kinase 15; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 7 protein; Cell division protein kinase 15; Serine/threonine-protein kinase ALS2CR7; Serine/threonine-protein kinase PFTAIR
Human Swiss Prot No.:	Q96Q40
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ALS2CR7 Polyclonal Antibody detects endogenous levels of ALS2CR7 protein.
Storage and Stability:	-20°C/1 year