

Antibodies Product Datasheet

Product Name:	CIB3 Antibody
Catalog Number:	C11530
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CIB3. AA range:138-187
Concentration:	1 mg/ml
Gene Name:	CIB3
Protein Name:	Calcium and integrin-binding family member 3
Other Name:	CIB3; KIP3; Calcium and integrin-binding family member 3; Kinase-interacting protein 3; KIP 3
Human Swiss Prot No.:	Q96Q77
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CIB3 Polyclonal Antibody detects endogenous levels of CIB3 protein.
Storage and Stability:	-20°C/1 year