

Antibodies Product Datasheet

Product Name:	DLC-1 Antibody
Catalog Number:	C11532
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RHG07. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	DLC1
Protein Name:	Rho GTPase-activating protein 7
Other Name:	DLC1; ARHGAP7; KIAA1723; STARD12; Rho GTPase-activating protein 7; Deleted in liver cancer 1 protein; DLC-1; HP protein; Rho-type GTPase-activating protein 7; START domain-containing protein 12; StARD12; StAR-related lipid transfer protein
Human Swiss Prot No.:	Q96QB1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DLC-1 Polyclonal Antibody detects endogenous levels of DLC-1 protein.
Storage and Stability:	-20°C/1 year