

Antibodies Product Datasheet

Product Name:	Mash2 Antibody
Catalog Number:	C11584
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Mash2 . at AA range: 1-80
Concentration:	1 mg/ml
Gene Name:	ASCL2
Protein Name:	Achaete-scute homolog 2
Other Name:	ASCL2; BHLHA45; HASH2; Achaete-scute homolog 2; ASH-2; hASH2; Class A basic helix-loop-helix protein 45; bHLHa45; Mash2
Human Swiss Prot No.:	Q99929
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Mash2 Polyclonal Antibody detects endogenous levels of Mash2 protein.
Storage and Stability:	-20°C/1 year